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**Stand der Information:** 06/2013 - V3.0

**Medizinprodukt:** Klasse II b, CE 0408

It is generally assumed that different types of wounds are colonized by microorganisms. This does not cause disease, but decontamination is always desirable in order to prevent wound infections and to provide optimal support for the wound healing process by cleaning the wound thoroughly and by sparing of viable tissue.

## General information – ActiMaris® Wound gel (20g)

for the decontamination, cleaning (biofilm) and moistening of acute and chronic wounds, skin, mucous membranes and first and second degree burns. It is also used as wound filler.

### 1. Product profile

ActiMaris® – Wound gel consists of sea salt, water and reactive oxygen.

Reactive oxygen eliminates wound typical pathogens, without damaging the tissue. The decontaminating effect against MRSA / VRE has been tested. The oxidative effect of singlet oxygen ( $O_2$ ) is also responsible for the neutralization of wound odours.

The alkaline pH - value of ActiMaris® – Wound gel additionally supports the controlled release of reactive oxygen, by providing the OH groups (as a part of the naturally balanced redox system).

Sea salt contains sodium chloride as its main ingredient and over 80 of the other different components of sea water. In addition to potassium, calcium, and magnesium, sulphur, bromine, iron and essential traces of iodine, zinc, manganese and selenium are also present. As in seawater, these components are present in their natural complex compounds and in physiologically correct concentrations so they favourably influence the cell milieu and (wound) metabolism.

The increased salinity of ActiMaris® - Wound gel allows for the optimal hyperosmotic effect so that swellings subside rapidly.

### 2. Fields of application

ActiMaris® - Wound gel is suitable:

#### For cleaning, reduction of swelling, moistening and decontamination of

- acute mechanical wounds (skin lacerations, bites, cuts, abrasions, tears and bruises) and postoperative wounds.
- chronic wounds (decubitus ulcers, arterial-venous ulcers, diabetic ulcers)
- thermal and chemical wounds (first and second degree burns)
- entry portals of urological catheters and PEG tubes and drains.

ActiMaris® - Wound gel is particularly suitable for cleaning of severe fibrinous or purulent, infected or strong foal-smelling wounds and as a filler of deep defect wounds.

### 3. Instructions for use

ActiMaris® - Wound gel should be applied at room temperature; if necessary, the wound gel can be warmed to body temperature before use. Before application of the wound gel, the wound should be cleaned with ActiMaris® - Wound Irrigation Solution (Forte or Sensitiv). The wound gel will be then generously applied to or in the wound and covered with an inert compress, gauze or similar. The duration and frequency of application should be adapted to the wound conditions and continued until all coatings and necrosis are easily removable and the wound is optically clean.

Depending on the wound condition, a bandage change can be done according to indications and phases, every 2 to 3 days. In case of non-infected wounds, necrosis and coatings can be mechanically removed depending on the wound condition.

In cases where bandaging is difficult to release, it is recommended to soak the wound bandages with ActiMaris® Wound rinsing solution (Forte or Sensitiv) for at least 5 minutes, so they can be gently removed without traumatizing the wound surface. In case of thick, adherent fibrinous or necrotic wound coatings, a primary surgical sanitation is recommended.

For the cleaning of the wound ActiMaris® - Wound gel can be used in combination with inert bandages or wound dressings (see also section 8: Interactions).

### 4. Tissue tolerance and biocompatibility

ActiMaris – Wound gel is suitable for long-term use even on sensitive / irritated skin and mucosa due its non-irritating and well-tolerated application. There is no inhibition of granulation and epithelialisation. Dermatologically approved.

### 5. Side effects

No side effects known at present.

### 6. Contraindications

ActiMaris® – wound gel is not to be applied in case of known or suspected allergy to any of the ingredients

### 7. Restrictions on use:

*Pregnancy and lactation:*

There are no indications of mutagenic or embryo toxic effects of the ingredients. Due to the lack of studies and long-term clinical experience with pregnant and breast feeding women, ActiMaris® - Wound gel should be used in these cases only after an accurate medical assessment.

*Infants and newborns:*

Due to the lack of studies and long-term clinical experience with pregnant and breast feeding women., ActiMaris® – Wound gel should be used in these cases only after an accurate medical assessment.

### 8. Interactions:

ActiMaris® - Wound gel works due to the oxidation of the reactive oxygen's content.

To avoid possible interactions, wound gel should not be mixed with or used in combination with substances containing silver, biguanides, taurolidine or iodine-containing solutions and wound dressings. This shall also apply even to wetting of bandages and wound dressings; here can only sterile, not active bandaging materials such as (e.g.) active agent free dressings, gauze, pads, sponges, gels, hydro fibres, alginates, hydrocolloids and similar should be used. Simultaneous application of wound gel with other ActiMaris® - products (Sensitiv and Forte – Wound rinsing solution) is possible.

### 9. General safety instructions

ActiMaris® – Wound gel is suitable only for external wounds on skin or mucosa. Do not use for infusion or injection! Do not swallow!

Use only closed and undamaged containers.

Keep out of reach of children!

It is possible for coloured fabrics to fade after contact with ActiMaris® - products.

### 10. Composition and technical information

Sea salt, water, lithium magnesium sodium silicate

Clear, colourless, basic and hypertonic gel with characteristic smell; contains 3% NaCl (Sal Maris) and 0,2% NaOCl.

### 11. Storage specification, duration, expiry date

ActiMaris® - Wound gel has to be stored at room temperature between +5° and +25°C.

Protect from direct sunlight.

Shelf life is 30 months. Expiry date is on the packaging.

The period of use is 3 months after opening the tube. Close carefully after use.

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ActiMaris AG **Fabbricante:**
ActiMaris AG